



Hohenstein Quality Label "Fit"

Objective:

With the help of the independent product endorsement "Best suited to the material and functional fit", you can guarantee a consistent fit for your customers – right across all your collections.

Description:

- The Hohenstein quality label is particularly suitable for NOS goods
- The set of sizes will be picked up on-site by our independent inspectors
- The fit is tested both in new condition and after washing as described on the care label
- The assessment is made by experienced clothing engineers using fitting models whose body measurements correspond to those in the "Size Germany" survey
- If the set of sizes meets all requirements, the product can be endorsed with the Hohenstein quality label at the POS

Requirements for test samples:

A product that reaches the highest assessment level both before and after washing can be awarded the Hohenstein quality label. This assures that the product satisfies the following requirements:

- The garment must fit with regard to all length and width measurements
- It must allow freedom of movement as appropriate for its purpose and function
- The fit and appearance must not be impaired by defects in workmanship
- The labels showing the used materials and care instructions must be present and correct in accordance with the Textile Labelling Act (EU) No. 1007/2011

The scope of service:

- Label printing template valid for 1 year
- Test report incl. fitting score
- If required, a clarifying meeting with try-on session

Customer benefit:

- The neutral assurance of product quality by an independent institution can be shown to the end user at the POS in the form of a hang-tag
- Quality assessment in new condition and after washing
- Effective prevention of complaints at the POS

What we need from you:

- A set of sizes in the original material incl. care instructions
- To extend the label for another year, our inspector picks up the latest set of sizes for re-testing before the expiry of the period of validity

Time required:

- Maximum 10 working days; express service: 5 working days

HOHENSTEIN INSTITUTE

Schloss Hohenstein
74357 Boennigheim

Contact

Clothing Technology Department
Stephanie Brenner
Head of Fitting
Tel. no.: +49 7143 271 304
Fax: +49 7143 271 94304
Email: s.brenner@hohenstein.de

www.hohenstein.de